INSTRUCTIONS FOR USE
SUCTION LINER SYSTEM

PRODUCT NAME
Suction Liner System

MODEL AND SPECIFICATION
YDI, YDII, YDIII, YDIV, YDV, YDVI

SHELF LIFE
3 years.

PRODUCT DESPCRIPTION
Suction liner system consists of trolley, canisters, negative pressure gauge and hose. It is intended for use in conjunction with suction liner, negative-pressure pump or central suction system. Re-use.

INTENDED PURPOSE
Suction liner system is intended for use in hospital clinical departments and during or after surgery to collect body fluids (blood, gastric fluid) and secretions (sputum, irrigation fluid) from patients.

INSTRUCTIONS
1. Before using, check the device for damage.
2. Insert the suction liners with required specification into the canisters and fix them, connect the suction liners in series in turn, and finally connect suction liner with the appropriate connection of the canister.
3. Connect the hose to the appropriate connection of the negative pressure pump or central suction system.
4. Turn on the negative pressure pump or central suction system, adjust the negative pressure to the desired pressure to use.
5. After use, keep the negative pressure, connect the free end of the soft tube on suction liner with the PATIENT end of the suction liner to make it airtight, turn off the negative pressure pump or central suction system, remove the suction liner. Do not carry the soft tube of suction liner.

PRECAUTIONS
1. DO NOT use the damaged device.
2. DO NOT put this device in direct contact with human wounds.

STORAGE AND TRANSPORTATION
1. Store suction liner system in a dry and ventilated room with relative humidity not higher than 80%, no corrosive gas.
2. During transportation, the product should be protected from sunlight, rain and snow, and must not be mixed with toxic or harmful products for transportation.

SYMBOLS ON LABEL

Symbol	Title of symbol	Symbol	Title of symbol	Symbol	Title of symbol
	Manufacturer		Use-by date		Keep dry
	Batch code		Non-sterile		Stacking limit
	Date of manufacture		CE Mark		Keep away from sunlight
	Fragile, handle with care		Consult instructions for use		Medical device
	Authorized representative in the European Union		Unique device identifier		Distributor
	Importer	-	-	-	-


Shijiazhuang Yuanding Medical Device Co., LTD
Building 102, No.2 Xingye Street, Economic Development Area, Xinle, Shijiazhuang City, 050700 Hebei, P.R. China
Tel/Fax: +86-311-88519708


Shanghai International Holding Corp. GmbH (Europe)
Eiffestrasse 80, 20537 Hamburg, Germany
Tel: +49-40-2513175 Fax: +49-40-255726